

# Motor API

Revision 1 - Software v2023.2.2

irw

## Summary

The following includes a description of the available methods used to interact with the motor controller via Python.

---

## Contents

|   |          |
|---|----------|
| <b>Arduino class</b>  | <b>2</b> |
| connect(self, port = None, baudrate = 19200) . . . . .          | 2        |
| close(self) . . . . .   | 2        |
| format(self, fragments) . . . . .                               | 2        |
| query(self, message) . . . . .                                  | 2        |
| check_connection(self) . . . . .                                | 2        |
| check_parameters(self) . . . . .                                | 3        |
| print_parameters(self) . . . . .                                | 3        |
| check_lobe_delays(self) . . . . .                               | 3        |
| print_lobe_delays(self) . . . . .                               | 3        |
| set_breath_count(self, n) . . . . .                             | 3        |
| set_delay(self, maneuver, delay) . . . . .                      | 3        |
| set_maneuver_order(self, maneuver) . . . . .                    | 4        |
| set_lobe_default(self, parameter, value, motors) . . . . .      | 4        |
| update_lobe_delays(self, value = None, motors = None) . . . . . | 4        |
| prepare_maneuver(self, maneuver, steps, motors) . . . . .       | 4        |
| run_maneuver(self) . . . . .                                    | 4        |
| run_profile_constant(self) . . . . .                            | 4        |
| run_profile_variable(self) . . . . .                            | 5        |

## Arduino class

Instantiate a controller of the `Arduino` class with the following:

```
controller = Arduino(port = 'COM1') # Windows  
controller = Arduino(port = '/dev/ttyACM0') # Linux
```

The available methods are listed below, including an example and description of each.

Note:

- Lobes are always specified in order of RU, RM, RL, LU, LL.
- Maneuvers may be specified as 'profile', 'inhale', or 'exhale'

---

### **connect(self, port = None, baudrate = 19200)**

```
controller.connect()
```

Connect to the Arduino at the specified port. If no port is specified, attempt to connect on the default port supplied with creating the `Arduino` object.

If successful, sets the `dev` instance variable from the serial connection.

---

### **close(self)**

```
controller.close()
```

Closes the connection to the motor controller.

---

### **format(self, fragments)**

```
controller.format(['SLS', '200', '11100'])
```

Formats a series of parameters to be read by the motor command protocol. The `fragments` argument can be a single string (e.g. `controller.format('SOI')`) or a list (e.g. `controller.format(['SDE', '5'])`).

Returns the formatted message but does not send to the controller.

---

### **query(self, message)**

```
controller.query('?S')
```

Sends a message to the motor controller.

Returns the message response from the motor controller.

---

### **check\_connection(self)**

```
controller.check_connection()
```

Tests the serial connection to the motor controller.

Returns the controller response.

---

### **check\_parameters(self)**

`controller.check_parameters()`

Query the motor controller global parameters (breath count, maneuver delays, default values).

Returns the controller response.

---

### **print\_parameters(self)**

`controller.print_parameters()`

Aesthetic printing of default controller parameters.

---

### **check\_lobe\_delays(self)**

`controller.check_lobe_delays()`

Query the delays stored in the current maneuver and in memory for each lobe/motor.

Returns the controller response.

---

### **print\_lobe\_delays(self)**

`controller.print_lobe_delays()`

Aesthetic printing of default controller parameters and delays.

---

### **set\_breath\_count(self, n)**

`controller.set_breath_count(5)`

Set the number of breaths to be run in a profile.

---

### **set\_delay(self, maneuver, delay)**

`controller.set_delay('inhale', 0.1)`

Set the delay for one of the global parameters. Maneuver options include:

- 'profile'
- 'inhale'
- 'exhale'

Returns the controller response.

---

**set\_maneuver\_order(self, maneuver)**

`controller.set_maneuver_order('inhale')`

Set the first maneuver of a breathing profile. Maneuver options include:

- 'inhale'
- 'exhale'

Returns the controller response.

---

**set\_lobe\_default(self, parameter, value, motors)**

`controller.set_lobe_default('delay', 1000, '11111')`

Set the default value for lobe/motor parameters. Parameter can be one of 'steps' or 'delay'. Step values can be 0-700, and delay values can be a positive integer less than 65535. Delay values are related to processor cycles, not a strict time value.

---

**update\_lobe\_delays(self, value = None, motors = None)**

`controller.update_lobe_delays()`

Sets the delays for each lobe/motor to the default value for all steps.

*Future method updates can interpret parameters in relation to the SAC command.*

---

**prepare\_maneuver(self, maneuver, steps, motors)**

`controller.prepare_maneuver('inhale', 200, '10111')`

Prepares an inhalation or exhalation maneuver for the defined step value and given motors.

Returns the controller response.

---

**run\_maneuver(self)**

`controller.run_maneuver()`

Runs the last prepared maneuver for each lobe.

Returns the controller response.

---

**run\_profile\_constant(self)**

`controller.run_profile_constant()`

Runs a breathing profile based on the global settings and constant-delay breathing maneuvers.

Returns the controller response.

---

**run\_profile\_variable(self)**

`controller.run_profile_variable()`

Runs a breathing profile based on the global settings and variable-delay breathing maneuvers, where the variable delay values are retrieved from memory.

---