

# Motor Command Protocol

Revision 2 - Software v2023.2.2

irw

## Summary

The following includes a description of the available commands used to control the motors.

---

## Contents

<b>Usage</b>	<b>2</b>
Main commands . . . . .	2
<b>General commands</b>	<b>3</b>
? . . . . .	3
?S . . . . .	3
?A . . . . .	3
<b>Global settings</b>	<b>4</b>
SN/X . . . . .	4
SD{IEP}/X . . . . .	4
SO{IE} . . . . .	4
<b>Lobe settings</b>	<b>5</b>
SLS/X/00000 . . . . .	5
SLD/X/00000 . . . . .	5
SM{IE}/X/00000 . . . . .	5
SAC . . . . .	5
SAC/X/00000 . . . . .	6
Scripting commands . . . . .	6
<b>Run commands</b>	<b>7</b>
RUN . . . . .	7
PROFILEC . . . . .	7
PROFILEV . . . . .	7

## Usage

- All commands are case-sensitive
- All commands must be enclosed in angle brackets (i.e. <Command>)
- Maximum command length is 32 characters between brackets
- The separator to be used within commands is /
- If a command has multiple parameter options, they are denoted within braces ({}). For example, the command SD{IEP}/X can take the form <SDI/X> or <SDP/X> to set respective values
- If a command is not defined, **Invalid command** will be returned
- A maximum of 700 consecutive steps can be executed per maneuver
- One motor corresponds to one lobe, and they are listed in a five-bit sequence (00000)

Lobe/Motor	Abbreviation	Bit
Right upper	(RU)	1----
Right middle	(RM)	-1---
Right lower	(RL)	--1--
Left upper	(LU)	---1-
Left lower	(LL)	----1

## Main commands

The main commands to use in a ready-to-run setting are for checking motor parameters, preparing a single maneuver, and running a profile (series of inhale/exhale maneuvers).

### Check commands

- ?
  - Query connection
- ?S
  - Query parameters

### Run commands

- SM{Maneuver}/{Steps}/{Motors}
  - Prepare a maneuver
- RUN
  - Run the prepared maneuver
- PROFILEC
  - Run a constant-delay breathing profile series. This uses a single delay value for both inhale and exhale maneuvers, set by the default value and number of steps for each lobe.
- PROFILEV
  - Run a variable-delay breathing profile series. This uses the specific values stored in memory for inhale and exhale delays. The number of steps is based on the default value.

## General commands

**?**

Query the connection to the controller.

### Input

None

### Output

Returns **OK-Motors** if the connection is working.

---

**?S**

Query the global and lobe-specific default settings.

### Input

None

### Output

Returns the following:

- Profile delay
  - Profile cycles
  - Motor default number of steps
  - Motor default delay
  - Motor step position (recorded since time on)
  - Inhale delay (before inhalation maneuver)
  - Exhale delay (before exhalation maneuver)
  - Starting maneuver (first maneuver in cycle: 1 for inhale, -1 for exhale)
- 

**?A**

Query all delays for each lobe.

### Input

None

### Output

Returns the following:

- Current delay values for each lobe
  - Inhale delay values stored in memory
  - Exhale delay values stored in memory
-

## Global settings

Format: `S{Setting}/{Value}/{Additional parameters}`

### SN/X

Set number. Sets the number of breathing cycles performed in a profile, a set of paired inhalation and exhalation maneuvers.

#### Inputs

- **X** - The number of cycles. Positive integer.

#### Outputs

None

---

### SD{IEP}/X

Set delay inhale/exhale/profile. Set the delay value taken before the inhalation maneuver, exhalation maneuver, or breathing profile, in seconds.

Takes one of the following forms: **SDI**, **SDE**, **SDP**.

#### Input

- **SD{IEP}**
  - **I** - inhalation
  - **E** - exhalation
  - **P** - profile (repeated breathing cycles)
- **X** - Number of seconds to pause before a maneuver
  - Float (decimal value) acceptable

#### Output

None

---

### SO{IE}

Set order inhale/exhale. Set the order of the maneuvers for profile breathing.

Takes one of the following forms: **SOI**, **SOE**.

#### Input

- **SO{IE}**
  - **I** - start with inhalation (1)
  - **E** - start with exhalation (-1)

#### Output

None

---

## Lobe settings

### SLS/X/00000

Set lobe steps. Sets the default number of steps to be taken by a given set of motors. Multiple motors can be set simultaneously.

#### Input

- X - Number of steps. Integer value. Max = 700.
- 00000 - Motor bits

#### Output

Prints “Not set” if a motor is not set by this command.

---

### SLD/X/00000

Set lobe delay. Set the default delay value for each in a given set of motors. Multiple motors can be set simultaneously.

#### Input

- X - Integer value of the delay. Arbitrary units related to processor cycles. (0-65535)
- 00000 - Motor bits

#### Output

Prints “Not set” if a motor is not set by this command.

---

### SM{IE}/X/00000

Set maneuver inhale/exhale. For a given set of motors, prepares a maneuver of X steps, in the given direction. Sending this command does not run the motion. This command should be followed by RUN to start a maneuver.

Takes one of the following forms: SMI/X/00000, SME/X/00000.

Sending an updated command with the same motor bit(s) set, before sending the RUN command, will invalidate previous maneuver settings for the same motor (i.e. only one maneuver can be set for a given motor at a time).

#### Input

- SM{IE} - Set the maneuver type (I - inhale/E - exhale)
  - X - Set the number of steps for the maneuver. Max = 700.
  - 00000 - Motor bits
- 

### SAC

Set all constant. Set the delay values for all motors and all steps to the default value for each motor.

### Input

None

### Output

None

---

### SAC/X/00000

A special case of the **SAC** command with parameters. For the given motors, set the supplied value as the delay value for all steps.

### Input

- **X** - Delay value (0-65535)

### Output

None

---

## Scripting commands

### SAI/I/X/00000

Set at index. For a given set of motors, set the delay value at an individual step. Check that the value was written correctly.

### Input

- **I** - index of the step to write
- **X** - value of the delay
- 00000 - motor bits

**Output** **OK** if the value was written correctly. **Error** if the value stored in the array does not match the value supplied.

---

### C/S/X/00000

*Experimental. Not for general use.*

Constant. Manually set constant value for delays. Set the supplied value as the delay for the given motors, up to the given number of steps.

### Input

- **S** - Number of steps. Integer value. Max = 700.
- **X** - Delay value (0-65535)
- 00000 - Motor bits

**Output** **OK** if function completes.

---

## Run commands

### **RUN**

Runs the single most recent maneuver stored for each motor via the **SM** command. Only runs maneuvers since the last **RUN** command. If multiple **SM** command have been sent, only the most recent is performed.

#### **Input**

None

#### **Output**

None

### **PROFILEC**

Runs a constant-delay breathing profile based on the global breathing parameters and default delay arrays for each motor.

Executes the following:

- Set constant-value delays
- Delay profile
- For the given number of breaths:
  - Prepare first maneuver (inhale/exhale)
  - Delay first maneuver (inhale/exhale)
  - Run maneuver and wait until completion
  - Prepare second maneuver (inhale/exhale)
  - Delay second maneuver (inhale/exhale)
  - Run maneuver and wait until completion

#### **Input**

None

#### **Output**

None

### **PROFILEV**

Runs a variable-delay breathing profile based on the global breathing parameters and Bezier-defined delay arrays for each motor.

Executes the following:

- Delay profile
- For the given number of breaths:
  - Prepare first maneuver variable-value delays (inhale/exhale, assigned to each lobe)
  - Prepare first maneuver (inhale/exhale)
  - Run maneuver and wait until completion
  - Prepare second maneuver variable-value delays (inhale/exhale, assigned to each lobe)
  - Prepare second maneuver (inhale/exhale)
  - Run maneuver and wait until completion

#### **Input**

None

## Output

None